


Proximate analysis of a coal

INTRODUCTION

The proximate analysis, allows classifying a coal and gives information for its industrial use. Coal is particularly classified according to its volatiles content : anthracite (8 to 18 %), semi-bituminous (18 to 26 %), bituminous (26 to 40 %), semi-bituminous (40 to 50 %). The thermogravimetric test allows giving a rapid result for a coal. All the operations heating and gas switching are automatically achieved by the controller associated with the thermoanalyzer.



EXPERIMENT

- Samples : coal
- Mass : 46.85 mg
- Crucible : Silica
- Atmosphere : Argon, then oxygen
- Heating mode :
Heating to 110°C and isothermal level during 5 minutes under argon
Rapid heating to 950°C and isothermal level during 10 minutes under argon
Switch from argon to oxygen at 950°C

RESULTS AND CONCLUSION

The proximate analysis of coal sample gives the following results :

- Moisture: 3.2 %
- Volatiles: 27.2 %
- Fixed carbon : 51.2 %
- Ash: 18.4 %

According to the volatiles content, the sample can be classified as a bituminous coal.

INSTRUMENT

THEMYS ONE TGA



- HIGH SENSITIVITY BALANCE FOR THE DETECTION OF SMALL MASS VARIATIONS specifically designed for TGA analysis.
- CONVENIENCE OF ONE FURNACE to reach temperatures as high as **1150°C or 1600°C**.
- PLUG AND PLAY INTERCHANGEABLE RODS to perform TGA only, TG-DSC, TG-DTA, and 3D high sensitivity/Cp measurements.
- EXTERNAL COUPLING CAPABILITY including evolved gas analysers