

Analysis of plaster with large volume TGA

INTRODUCTION


Large volume thermogravimetry may be of great interest, because it makes possible the analysis of non-homogeneous samples. It is especially interesting in the case of natural raw materials or ores, because with a larger sample size, the sampling becomes less critical.

Using large volumes is also a way to obtain more accurate measurements : it can be especially interesting if the mass variation is less than 1%.

With the THEMYS LV, it is possible to use crucibles with a volume of 18 mL. (diameter = 22 mm, height = 50 mm)

EXPERIMENT

- A sample of plaster was analyzed in TGA mode .
- Sample mass : 13.35 g
- The temperature was programmed from ambient to 300°C at 1 K/min.
- Atmosphere : air




RESULTS AND CONCLUSION

The figure shows that we can separate the mass loss in two steps : between ambient and 145°C : a mass loss of 0.373 % is measured : it is due to the dehydration of gypsum.

Between 145°C and 250°C : a mass loss of 5.763 % is measured. It is due to the dehydration of plaster.

INSTRUMENT

THEMYS LV TGA



- **EXTERNAL COUPLING CAPABILITY**
designed for evolved gas analyzers (FTIR, MS, GCMS, MSFTIR, or FTIR-GCMS)
- **ULTRA-HIGH TEMPERATURE CAPABILITY**
to 2000°C with a single furnace.
- **MODULAR ADAPTIONS ALLOWING**
TGA only, DTA only, TG-DTA up to 2000°C, High temperature calorimetry up to 1600°C all in one instrument
- **HIGH ACCURACY & VERSATILITY**
hang-down symmetrical beam balance, specifically designed for TGA applications
- **ACCURATE AND SENSITIVE ULTRA-HIGH TEMPERATURE**
heat flow measurement with Tri- Couple DTA technology