

# **PROCESS** SAFETY

# THERMAL ANALYSIS, CALORIMETRY & THERMOKINETICS SOLUTIONS

QUICK STABILITY SCREENING •
THERMAL STABILITY WITH PRESSURE AND EVOLVED GAS DATA •
PROCESS / SYNTHESIS REACTION UNDERSTANDING •
ADIABATIC / ACCELERATING RATE CALORIMETRY •
MODELLING AND PREDICTION WITH KINETICS ANALYSIS SOFTWARE •



In chemical, pharmaceutical and many more industrial settings process safety is of critical importance. Thermal risks must be assessed in normal and runaway process conditions to avoid costly downtime and damaging incidents to health, the environment and corporate reputations.

The Process Safety laboratory is fundamental to any chemical facility. Whether synthesising small or large amounts of material each process step and compound needs assessing for thermal stability.

To avoid thermal hazards the entire process from lab, to pilot, to plant, needs safe scale-up support with dedicated instruments. KEP Technologies understands your challenges and offers a choice of solutions that provide experimental control, instrument versatility and quality results.

### COMMON PROCESS SAFETY STUDIES & SOLUTIONS



DSC stage and due to the central role of gas releasing reactions in plant process hazards, larger samples of about 5-10g are heated to destruction in cells equipped with a pressure sensor. The focus of the study is to understand volume of gas evolved, when and how fast it's evolved in order to understand it's decomposition profile.

### Thermal Stability with Pressure and Evolved Gas

Adiabatic conditions are a worst case scenario that assumes no heat loss from the process, but are quite realistic for large plant vessels that lose heat slowly. Here an adiabatic reaction calorimeter (ARC) tests samples at minor to no heat loss conditions, and measures the largest temperature rise achieved to establish the worst case scenario if decomposition occurs. Depending on process specifics, a low phi adiabatic study can also be done to more closely mimic large plant conditions and determine vent sizing.

Adiabatic Tests

#### Using small samples a fast indication of thermal instability can be gained by running all process materials through a DSC Thermal Analyzer. The DSC measures the temperature and heat of decomposition, with the main study focus being when and how much heat is given off. As sample sizes are small and pressure data isn't an output of DSC screening, if potential process safety issues are indicated, such as exotherms, further study is required.

Quick Screening

Once decomposition risks and boundaries are established, it's necessary to understand the desired chemical reaction. Reaction Calorimetry determines how much heat is given out under normal operating conditions. The focus here is on heats of reaction, heat flows and gas evolution, with study results being the heat of reaction and heat release rate. The cooling capacity of the reactor should also be checked to be sufficient to control the reaction temperature profile from the prior step and so avoid runaway situations.

### **Process Understanding**

The best software serve as a complement to all of the above studies by modelling scenarios to reduce the number of actual experiments, vastly

speeding-up process safety testing.

Modelling with software

١

١

T

"If I was setting up a new process safety laboratory, I would start with a CALVET Calorimeter because even if I couldn't solve each process safety problem, it can be used for both synthesis reactions (mixing cell) and decomposition reactions, by offering an appreciable advantage, namely pressure measurement. I would start with a CALVET even if it meant supplementing it later on."

Professor Stoessel - TÜV SÜD Process Safety Basel, Switzerland

### THE KEP TECHNOLOGIES ADVANTAGE

KEP Technologies is radically step-changing it's coverage of the process safety market by offering the widest and most versatile choice of solutions. Now you can consult with one company, KEP Technologies, to address your challenges across the broadest number of process safety studies on the market.

Each solution embodies our "Reimagine Material Characterization" value proposition by delivering strongly against the three core customer benefits of **Experimental Control**, **Instrument Versatilty and Quality Results**.

We believe solutions that provide these benefits will deliver the highest value to our customers.

In addition to our core range we are able to provide **customized solutions** by harnessing the engineering and project management of our highly skilled organisation.



### **CUSTOMIZED SOLUTIONS**

Modular design allowing for upgraded and tailored functionality Access to all previous non-proprietary custom requests Open access to engineering development team



### **INSTRUMENT**

## SETLINE



### **SPECIFICATIONS**

Temperature range (°C)	-170 to 700
Programmable heating rate (°C/min)	0.01 to 100
Enthalpy accuracy / precision* (%)	+/- 0.8 / 2.5
Temperature accuracy / precision* (°C)	+/- 0.3 / 0.50
DSC measurement range (mW)	+/- 6 000
Atmosphere	Inert gas, air, High pressure crucibles up to 500 bar at 600 °C
Autosampler	SETLINE DSC+ version featuring a 59 positions autosampler

\* Based on indium melting tests

For more information on specifications please consult the product information and brochures available on our website : www.setaramsolutions.com

### SCREENING DSC : ANALYSIS OF PROPERGOL

### **INTRODUCTION**

The thermal decomposition of Propergol was studied in a tightly closed high pressure crucible and in an open crucible with milligram-scale samples.

#### **EXPERIMENT**

Both experiments were run at a scanning rate of 3°C/min, under inert gas conditions.

### **RESULTS AND CONCLUSION**

It was observed that if the pressure rises during decomposition (closed crucibles), both total heat and reaction schemes are different. However, three peaks remain at similar positions and shape (cf. arrows), which was confirmed by built-in deconvolution module of Calisto data treatment software.

The flexibility of Setline DSC gives fast insight into the behavior of chemicals under varying gaseous conditions (ex: ambient pressure vs. high pressure, oxidizing vs. inert conditions, etc...).



# THERMAL STABILITY WITH PRESSURE AND EVOLVED GAS

### **INSTRUMENT**



### **SPECIFICATIONS**

Temperature range (°C)	Room Temperature to 400
Temperature accuracy / control (°C)	0.01
Heating rate range (°C/min)	0.5 to 10
Modes	Temperature Scanning, Isothermal, Dual Scan
Pressure Range (bar)	0 to 200
Pressure Resolution (bar)	0.001
Pressure Accuracy (bar)	+/- 2

For more information on specifications please consult the product information and brochures available on our website : www.setaramsolutions.com

### **APPLICATION**

### **DECOMPOSITION OF DTBP BY RAPID SCREENING CALORIMETRY**

### **INTRODUCTION**

Di-Ter Butyl Peroxide (DTBP) is an unstable chemical used in the polymer industry to initiate polymerization reactions. As large quantities of such a product need to be stored at the plant, it is necessary to assess the risks with the thermal decomposition of this chemical. DTBP is also a typical compound used to assess the performance of Accelerating Rate Calorimeters and Rapid Screening Calorimeters.

### **EXPERIMENT**

The temperature and pressure increase data obtained are very repeatable.

### **RESULTS AND CONCLUSION**

5 g samples of the same 15 wt% (DTBP) solution in toluene were heated in a 8 mL titanium cell at 2 °C/min from 40 °C up to 300 °C.



# QUICK SCREENING THERMAL STABILITY WITH PRESSURE AND EVOLVED GAS PROCESS UNDERSTANDING

### **INSTRUMENT**



### **CALVET Calorimeter**

# **ISOTHERMAL OR TEMPERATURE SCANNING**

with Calvet 3D sensor capturing 93-95% of all heat forms, the highest level on the market

### **CONVENIENT, INTERCHANGEABLE CELLS to**

perform varied experiments with one instrument :

- Includes but is not limited to: mixing, high pressure and interactions between solids, liquids and gases
- Pressure measurement allows it to be used for both synthesis and decomposition reactions

### FOR SMALL SCALE REACTION CALORIMETRY

negating the need to use the larger reaction

with low temperature version CALVET CRYO and high temperature version CALVET HT

	CALVET	CALVET CRYO	CALVET HT	
Temperature range (°C)	Ambient to 300	-196 to 200	Ambient to 600	
Temperature accuracy (°C)	+/-0.3 *	+/-0.5 **	+/-1*	
Temperature precision (°C)	+/-0.15*	+/-0.25**	+/-0.5*	
Programmable temperature scanning rate	0.001 to 2°C/min	0.01 to 1°C/min	0.01 to 2°C/min	
Enthalpy accuracy	+/-0.4 *	+/-0.2 **	+/-1*	
Calorimetric precision (%)	+/-0.4*	+/-0.5**	+/-1.5*	
Cells (ml)	Up to 12.5 (standard cell)	Up to 12.5 (standard cell)	Up to 7	
Pressure measured and controlled (bar [psi])	350 [5,075]; 600 [8,700]; 1000 [14,600]	100 [1,450]; 600 [8,700]; 1000 [14,600]	100 [1,450]; 300 [4,350]; 400 [5,800]	

### **SPECIFICATIONS**

\* Based on indium melting tests \*\* Based on naphthalene melting tests

For more information on specifications please consult the product information and brochures available on our website : www.setaramsolutions.com

### ANALYSIS OF DIETHYL SULFATE

### **INTRODUCTION**

The thermal decomposition of diethyl sulfate was studied in a tightly closed high pressure cell equipped with a pressure measurement system.

### **EXPERIMENT**

The experiment was run at a scanning rate of  $0.5^{\circ}$ C/min.



### **RESULTS AND CONCLUSION**

It was observed that the biggest contribution to the pressure increase was due to a first decomposition at about 175°C. The second decomposition peak is probably consuming part of the gas produced at the first stage. More information can be extracted from this signal like pressure release rate, condensable / non condensable gas ratio.

### **POLYMERIZATION STUDY**

### **INTRODUCTION**

The polymerization reaction of Vinyl Pyrrolidone in presence of 4, 4'-azobis-cyanovaleric acid was studied using a CALVET Calorimeter with membrane mixing cells.

### **EXPERIMENT**

Vinyl Pyrrolidone and 4, 4'-azobis-cyanovaleric acid were placed in the two chambers of the cell, separated by a thin membrane. The calorimeter was run under isothermal mode at 50 °C. In-situ mixing was provided by piercing the membrane.

### **RESULTS AND CONCLUSION**

0.06 ∧ Exa |-0.06 0.05 0.05 0.04 0.04 HeatFlow (Mig) Heat : -685.6 (J/g) (D/M) 0.03 Peak 0.02 0.02 Heat : -85.4 (J/g) 0.01 0.01 0 0 2 5 Time (h) 10 3 4 6

It was observed that a first, sharp peak, probably linked with the initiation of the reaction, was followed by a slower kinetics and higher heat process. This second peak is linked with the polymerization of Vinyl Pyrrolidone.

Deconvolution of these peaks gives the heat of the initiation, and by difference, the heat of polymerization could be calculated.



### **INSTRUMENT**



### **SPECIFICATIONS**

Temperature range (°C)	Room temperature to 500
Temperature increase detection Threshold (°C/min)	0.005 to 0.02
Temperature Precision (repeatability, °C)	+/- 0.05
Modes	Heat-Wait-Search, Isothermal, Temperature scanning
Pressure Range (bar)	0 to 200
Pressure Resolution (bar)	0.001
Pressure Accuracy (bar)	+/- 2

For more information on specifications please consult the product information and brochures available on our website : www.setaramsolutions.com

### DECOMPOSITION OF DTBP USING ACCELERATED RATE CALORIMETRY

### **INTRODUCTION**

Peroxides, including Di-Ter Butyl Peroxide (DTBP), are typically unstable chemicals that require careful safety studies before being involved in industrial processes.

### **EXPERIMENT**

The following were heated in 8 mL titanium cells using the Heat-Wait-Search mode:

- 5 g samples of the same 15 wt% DTBP solution in toluene
- One 5 g sample of a 40 w% DTBP solution in toluene

The Heat-Wait Search parameters were:

- Start temperature: 97 °C
- Temperature steps: 5°C
- Soak Time: 30 min, wait time: 30 min, search time: 15 min
- Detection threshold: 0.02 °C/min
- End Temperature: 250 °C (400 °C with the 40 wt% DTBP solution)

### **RESULTS AND CONCLUSION**

The analysis of experimental data allows for the determination of the onset temperature of decomposition, the adiabatic temperature rise (raw and phi-factor corrected), the temperatures at maximum temperature rise and pressure rise rates, and the pressure increase under adiabatic conditions.





This series of tests has shown the impact of concentration on the temperature and pressure rise of DTBP. A significant increase of concentration leads to drastically higher thermal and pressure risks.

Accelerating Rate Calorimetry provides the necessary data to evaluate these risks, and the precision (repeatability) of TAC-500 AS measurements has been proven.

	Tonset (°C)	ΔTad, raw (°C)	ΔTad, corrected (°C)	ΔT at max T rate (°C)	ΔT at max P rate (°C)	ΔTad (Mpa)
DTBP 15% -1	121.56	56.17	104.50	169.16	164.73	1.90
DTBP 15% -2	121.70	54.46	101.32	166.34	172.07	1.84
DTBP 15% -3	121.64	54.89	101.90	164.30	168.35	1.83
DTBP 40%	116.79	138.71	257.44	210.55	187.46	9.70

# **MODELLING WITH SOFTWARE**

### **SOFTWARE**

## **AKTS Thermal Safety Software**



**KEP TECHNOLOGIES ALSO PROVIDES AKTS SOFTWARE** for precise modelling of runaway situations

### WITH AKTS CALORIMETRY DATA SETS provides

 the Self-Accelerated Decomposition
Temperature (SADT), Time to Maximum Rate (TMR) and the full package of thermal hazards data can all be precisely predicted

ADIABATIC CALORIMETER can be used as single experiment verification instead of for a series of experiments at the start of process evaluation as such reducing countless time and resources

For more information please consult the product information and brochures available on our website : www.setaramsolutions.com

### **DECOMPOSITION OF 3-METHYL-4-NITROPHENOL USING DSC AND CALVET**

The decomposition of 3-methyl-4-nitrophenol was studied at different heating rates using DSC and CALVET (A).

The experiments at different heating rates were treated with AKTS software (B).

The variation of the runaway time under true adiabatic mode (Phi Factor = 1) can be calculated for any process temperature (C). The critical value TMRad = 24 hours is obtained at  $153^{\circ}$ C in that case. Dashed lines depict the confidence interval of the calculation.



An adiabatic experiment with a Phi factor = 3.2 was performed for the final validation of the simulation, and compared to calculated adiabatic data (D).

SADT can be determined applying «Finite Element Analysis» (D).



Switzerland – France – China – United States – India – Hong Kong For contact details: www.setaramsolutions.com or setaram@kep-technologies.com

Setaram is a registered trademark of KEP Technologies Group