

Heating and cooling curve up to 1200°C

INTRODUCTION

The thermal characterization of materials can require fast heating and cooling cycles, so that it replicates the demanding conditions encountered in their lifespan. THEMYS FLASH is a Thermogravimetric Analyzer capable of replicating fast temperature cycles on samples of various materials, up to 600 °C/min..

EXPERIMENT

THEMYS FLASH is composed of 5 cavities equipped with a high-efficiency image furnace for accurate and responsive heating and cooling using infrared radiations. It was programmed to :


- heat up from 240°C to 1200°C
- stay isothermal during 25 minutes
- cool down to 480°C

RESULTS AND CONCLUSION

It takes less than 2 minutes for the furnace to reach 1000°C and less than 5 minutes 30 seconds to hit the target of 1200°C.

During the 25 minutes isothermal sequence, the temperature is kept at the setpoint within +/- 0.15°C.

During the cooling phase, it takes less than 40 seconds for the furnace to reach 1000°C, and about 7 minutes to cool down to 480°C.



INSTRUMENT

THEMYS FLASH

Ambient to 1200°C

