

## Solid phase transition in Zirconia

### INTRODUCTION

Phase transitions are commonly characterized using DTA or DSC sensors. These techniques are very sensitive at low temperatures, but may face certain limitations above 1000°C. SETARAM's THEMYS thermal analyzers offer an interesting alternative in the form of TRI-COUPLES DTA sensors. These sensors are designed with **3 thermocouple junctions** under each crucible instead of just one for standard DTA. They also use less material than a DSC sensor and are therefore less subject to parasitic radiative effects at high temperatures.

As a result, this design considerably improves the amount of heat collected and the sensitivity of the sensor.

Low energy phenomena at high temperature can be accurately measured with this sensor. For example, **zirconia solid-solid phase transition (Monoclinic - Tetragonal) around 1200°C** has been tested and compared with standard DSC sensor.

### EXPERIMENT

The following experimental procedure has to be used:

- Sample: Zirconia powder
- Mass: 83mg (DTA) 87mg (DSC)
- Reference: Empty Al<sub>2</sub>O<sub>3</sub> crucible
- Heating: 800 – 1500°C at 10K.min<sup>-1</sup>
- Apparatus: THEMYS with Type S tri-couple DTA or DSC sensor
- Atmosphere: Air at 20ml.min<sup>-1</sup>




Figure 1 : DTA curve of zirconia sample between 850°C and 1450°C, measured with DSC and DTA Tri-couple sensor

### RESULTS AND CONCLUSION

Upon heating, an endothermic effect is measured between 1125°C and 1225°C. It corresponds to the solid-state transition from monoclinic to tetragonal. While it is observed with both DTA Tri-couple and DSC sensor, **the sensitivity is 23% higher with the DTA Tri-couple**. Therefore, this sensor design allows a better measurement of thermal phenomena at high temperature. Recrystallization to monoclinic structure is also measured upon cooling, at much lower temperature.

### INSTRUMENT

#### THEMYS STA



**ACCURATE AND SENSITIVE ULTRA-HIGH TEMPERATURE**  
Heat flow measurement with Tri-Couple DTA technology

**ULTRA-HIGH TEMPERATURE CAPABILITY**  
to 2400°C with a single furnace.

**MODULAR ADAPTIONS ALLOWING**  
TGA only, DTA only, TG-DTA, and TMA up to 2400°C,  
DSC only and TG-DSC up to 1600°C all in one  
instrument

**EXTERNAL COUPLING CAPABILITY**  
designed for evolved gas analyzers (FTIR, MS, GCMS,  
MSFTIR, or FTIR-GCMS)