

# FLEXI CHILL



## POWERFUL AIR COOLED CHILLER FOR MULTIPLE LABORATORY OR INDUSTRIAL APPLICATIONS

Finest design of heat exchangers and cooling fans for high cooling power capabilities

## PLUG&PLAY, SIMPLICITY, EASY SETTINGS

- Effortless use with no temperature settings
- Limited maintenance with no refrigerant circuit, i.e. no leaks and potential environmental problems
- Fast connection to various systems or instrumentation

## ROBUST DESIGN

Developed in Switzerland by our research and innovation team and CE marked

PERFORMANCE	
<b>Cooling capacity</b>	1800 W at RT = 20°C 1250 W at RT = 30°C
<b>Pump flowrate</b>	up to 4L/min
<b>Pump outlet pressure</b>	Up to 1.5 bar
<b>Tank maximum capacity</b>	5 L
<b>Temperature range</b>	RT to 70°C
GENERAL	
<b>Size (W x D x H)</b>	420 x 420 x 370 mm 16.5 x 16.5 x 14.6 inch
<b>Weight</b>	15 kg 33 lb
<b>Power supply</b>	110 / 230 V 50/60 Hz